




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale PBJ, 4,200 kg

€1,900.00
plus VAT &
Shipping

Product Images



Description

- pBJ only: Internal automatic calibration when temperature changes and time-controlled according to defined intervals, guarantees high accuracy and makes location independent
- pBS only: CAL adjustment program for setting the accuracy
- Capacity display: An increasing luminous band indicates the still available weighing range
- Metal housing: rugged and stable
- Dosing aid: High-stability mode and other filter settings selectable
- Weighing with tolerance range (checkweighing): Input of an upper/lower limit value. A visual signal assists portioning, dosing and sorting
- Totalizing of weight values
- Balance identification number: 4-digit, freely programmable, is printed in the calibration record
- Automatic data output to PC/printer after each balance stability
- Draft shield included as standard, for models with readability = 0.001 g, weighing chamber W×D×H 180×193×87 mm
- Dust cover included in delivery

- Single-cell advanced technology:
- Fully automatic production of the load cell from one piece
- Stable temperature behavior
- Short settling time: 3 sec
- High mechanical robustness
- High corner load safety
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model. This code is made up of two quality features and, in addition to the product features, technical data and pictograms, provides you with a further decision-making aid for the perfect scale for your application.

Additional Information

No.	KP-5313
Manufacturer (Brand)	KERN
EAN	4058072148447
Transport temperature	Room temperature
Color	light greylgrey
Readability	10 mg
Operation via	Membrane keypad
Weight	6.5 kg
Reproducibility	10 mg
Scales features	with switchable weighing units
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Precision scale
Weighing range MAX	4 kg

