



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN precision scale PBJ, 620 g

**€2,050.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Internal automatic calibration with temperature changes and time-controlled according to defined intervals, guarantees high accuracy and makes independent of location only PBS: CAL calibration program for setting the accuracy Capacity display: A rising luminous band indicates the still available weighing range Metal housing: Rugged and stable Dosing aid: High-stability mode and other filter settings selectable Weighing with tolerance range (checkweighing): Input of an upper/lower limit value. An optical signal supports portioning, dosing and sorting Totalizing weight values Balance identification number: 4-digit, freely programmable, printed in the calibration record Automatic data output to the PC/printer after each balance stability Draft shield as standard, for models with readability = 0.001 g, weighing chamber W×D×H 180×193×87 mm Dust cover included Single-cell advanced technology: Fully automatic production of the load cell from one piece Stable temperature behavior Short settling time: 3 sec High mechanical ruggedness High corner load safety

## Additional Information

---

No.	KP-5315
Manufacturer (Brand)	KERN
EAN	4058072148461
Transport temperature	Room temperature
Color	light greylgrey
Readability	1 mg
Operation via	Membrane keypad
Weight	5.82 kg
Reproducibility	1 mg
Scales features	with switchable weighing units
Material weighing plate	Stainless steel
Weighing range MAX	620 g
Manufacturer	CORE
Type scale	Precision scale

