




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale PEJ, 2,200 kg

**Price on
request
plus VAT &
Shipping**

Product Images



Description

- PEJ: Internal automatic calibration/adjustment, guarantees high accuracy and makes it independent of location
- PES: CAL adjustment program for setting the accuracy using an external test weight
- Metal housing: rugged and stable
- Weighing with tolerance range (checkweighing): Input of 2 upper and 2 lower limits via 4 arrow keys. An acoustic signal assists portioning, dosing or sorting
- Draft shield standard on models with readability = 0.001 g. Weighing chamber W×D×H 170×150×100 mm.
- Below-balance weighing: Possibility of load suspension on the bottom of the balance by means of built-in eyelet, standard on models with readability = 0.001 g, hook not included in delivery
- Protective cover included in delivery
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model. This code consists of two quality features and provides you with a further decision-making aid for the perfect scale for your application, in addition to the product features, technical data and pictograms.

Additional Information

No.	KP-5352
Manufacturer (Brand)	KERN
EAN	4058072148805
Transport temperature	Room temperature
Color	light grey grey
Scales features	with switchable weighing units
Readability	10 mg
Operation via	Membrane keypad
Weight	7.5 kg
Manufacturer	CORE
Material weighing plate	Stainless steel
Reproducibility	10 mg
Type scale	Precision scale
Weighing range MAX	2 kg

