



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale PES, 15 kg

**Price on
request
plus VAT &
Shipping**

Product Images



Description

- PEJ: Internal automatic calibration/adjustment, guarantees high accuracy and makes the system location-independent
- PES: CAL adjustment program for setting the accuracy using an external test weight
- Metal housing: rugged and stable
- Weighing with tolerance range (checkweighing): Input of 2 upper and 2 lower limits via 4 arrow keys. An acoustic signal assists portioning, dosing or sorting
- Draft shield standard on models with readability = 0.001 g. Weighing chamber W×D×H 170×150×100 mm.
- Below-balance weighing: Possibility of load suspension on the bottom of the balance by means of built-in eyelet, standard on models with readability = 0.001 g, hook not included in delivery
- Protective cover included in delivery
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model. This code is made up of two quality features and, in addition to the product features, technical data and pictograms, provides you with a further decision-making aid for the perfect scale for your application.

Additional Information

No.	KP-5355
Manufacturer (Brand)	KERN
EAN	4058072148836
Transport temperature	Room temperature
Color	light greylgrey
Readability	100 mg
Operation via	Membrane keypad
Weight	5.69 kg
Reproducibility	100 mg
Scales features	with switchable weighing units
Material weighing plate	Stainless steel
Weighing range MAX	15 kg
Manufacturer	CORE
Type scale	Precision scale

