



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale PES, 4,200 kg

**Price on
request
plus VAT &
Shipping**

Product Images



Description

- PEJ: Internal automatic calibration/adjustment, guarantees high accuracy and makes the system location-independent
- PES: CAL adjustment program for setting the accuracy using an external test weight
- Metal housing: rugged and stable
- Weighing with tolerance range (checkweighing): Input of 2 upper and 2 lower limits via 4 arrow keys. An acoustic signal assists portioning, dosing or sorting
- Draft shield standard on models with readability = 0.001 g. Weighing chamber W×D×H 170×150×100 mm.
- Below-balance weighing: Possibility of load suspension on the bottom of the balance by means of built-in eyelet, standard on models with readability = 0.001 g, hook not included in delivery
- Protective cover included in delivery
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model. This code consists of two quality features and provides you with a further decision-making aid for the perfect scale for your application, in addition to the product features, technical data and pictograms.

Additional Information

No.	KP-5358
Manufacturer (Brand)	KERN
EAN	4058072148867
Transport temperature	Room temperature
Color	light greylgrey
Readability	10 mg
Operation via	Membrane keypad
Weight	6 kg
Reproducibility	10 mg
Scales features	with switchable weighing units
Material weighing plate	Stainless steel
Weighing range MAX	4 kg
Manufacturer	CORE
Type scale	Precision scale

