



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN precision scale PFB, 120 g

**€370.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Convenient operating philosophy: all main functions have their own button on the control panel
- Compact dimensions, advantageous where space is limited
- Capacity display: A rising luminous band indicates the remaining available weighing capacity
- Level indicator and foot screws for precise leveling of the balance as standard, ensuring the most accurate weighing results
- Draft shield as standard (only for models with weighing plate size  $\varnothing$  80 mm and 120 mm). Weighing chamber W×D×H 158×143×64 mm
- Dust cover included in delivery
- To help you make the right choice for your KERN scale, we have listed a quality code for each model. This code is made up of two quality features and, in addition to the product features, technical data and pictograms, provides you with a further decision-making aid for the perfect scale for your application.

## Additional Information

---

|                         |                                |
|-------------------------|--------------------------------|
| No.                     | KP-5366                        |
| Manufacturer (Brand)    | KERN                           |
| VGKL number             | 105837123                      |
| EAN                     | 4058072148942                  |
| Transport temperature   | Room temperature               |
| Color                   | light greylgrey                |
| Readability             | 1 mg                           |
| Operation via           | Membrane keypad                |
| Weight                  | 2.32 kg                        |
| Reproducibility         | 1 mg                           |
| Scales features         | with switchable weighing units |
| Material weighing plate | Stainless steel                |
| Weighing range MAX      | 120 g                          |
| Manufacturer            | CORE                           |
| Type scale              | Precision scale                |

