




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale PFB, 200 g

€390.00
plus VAT &
Shipping

Product Images



Description

- Convenient operating philosophy: all main functions have their own button on the control panel
- Compact dimensions, advantageous where space is limited
- Capacity display: A rising luminous band indicates the remaining available weighing capacity
- Level indicator and foot screws for precise leveling of the balance as standard, ensuring the most accurate weighing results
- Draft shield as standard (only for models with weighing plate size \varnothing 80 mm and 120 mm). Weighing chamber W×D×H 158×143×64 mm
- Dust cover included in delivery
- To help you make the right choice for your KERN scale, we have listed a quality code for each model. This code is made up of two quality features and, in addition to the product features, technical data and pictograms, provides you with a further decision-making aid for the perfect scale for your application.

Additional Information

No.	KP-5368
Manufacturer (Brand)	KERN
VGKL number	105837203
EAN	4058072148966
Transport temperature	Room temperature
Color	light greylgrey
Readability	1 mg
Operation via	Membrane keypad
Weight	2.3 kg
Reproducibility	2 mg
Scales features	with switchable weighing units
Material weighing plate	Stainless steel
Weighing range MAX	200 g
Manufacturer	CORE
Type scale	Precision scale

