




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN precision scale PFB, 6 kg

**€450.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Convenient operating philosophy: all main functions have their own button on the control panel
- Compact dimensions, advantageous where space is limited
- Capacity display: A rising luminous band indicates the remaining available weighing capacity
- Level indicator and foot screws for precise leveling of the balance as standard, ensuring the most accurate weighing results
- Draft shield as standard (only for models with weighing plate size  $\varnothing$  80 mm and 120 mm). Weighing chamber W×D×H 158×143×64 mm
- Dust cover included in delivery
- To help you make the right choice for your KERN scale, we have listed a quality code for each model. This code is made up of two quality features and, in addition to the product features, technical data and pictograms, provides you with a further decision-making aid for the perfect scale for your application.

## Additional Information

---

No.	KP-5374
Manufacturer (Brand)	KERN
EAN	4058072149017
Transport temperature	Room temperature
Color	light grey grey
Readability	50 mg
Operation via	Membrane keypad
Weight	2.69 kg
Reproducibility	50 mg
Scales features	with switchable weighing units
Material weighing plate	Stainless steel
Weighing range MAX	6 kg
Manufacturer	CORE
Type scale	Precision scale

