




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale PLJ, 6,200 kg

€1,280.00
plus VAT &
Shipping

Product Images



Description

- PLJ 2000-3A high-quality milligram balance with enormous weighing capacity up to 2100 g - ideal for large samples or heavy tare vessels
- Ergonomically optimized control panel for left- and right-handers
- KERN PLJ: internal automatic adjustment, guarantees high accuracy and makes independent of location.
- Ideal for mobile, legal-for-trade applications, such as ambulatory gold and jewelry purchasing
- KERN PLS: CAL adjustment program for setting accuracy, external test weights at extra cost
- Draft shield standard on models with weighing plate size Ø 110 mm, weighing chamber Ø×H 60×150 mm
- KERN PLS/PLJ-F: Strain gauges
- KERN PLS/PLS-A: force compensation
- Protective cover included in delivery
- Plug-in power supply KERN YKA-05 suitable for PLS-A and PLJ-A
- Plug-in power supply KERN YKA-04 suitable for PLS-F and PLJ-F
- To make the right choice of your KERN scale easier for you, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect scale for your application.

Additional Information

No.	KP-5395
Manufacturer (Brand)	KERN
VGKL number	105834620
EAN	4058072149178
Transport temperature	Room temperature
Color	light greylgrey
Readability	10 mg
Operation via	Membrane keypad
Weight	7.04 kg
Reproducibility	10 mg
Scales features	with switchable weighing units
Material weighing plate	Stainless steel
Weighing range MAX	6 kg
Manufacturer	CORE
Type scale	Precision scale

