



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale PLJ, 6,200 kg

€1,330.00
plus VAT &
Shipping

Product Images



Description

- PLJ 2000-3A high-quality milligram balance with enormous weighing capacity up to 2100 g - ideal for large samples or heavy tare vessels
- Ergonomically optimized control panel for left- and right-handers
- KERN PLJ: internal automatic adjustment, guarantees high accuracy and makes independent of location.
- Ideal for mobile, legal-for-trade applications, such as ambulatory gold and jewelry purchasing
- KERN PLS: CAL adjustment program for setting accuracy, external test weights at extra cost
- Draft shield standard on models with weighing plate size Ø 110 mm, weighing chamber Ø×H 60×150 mm
- KERN PLS/PLJ-F: Strain gauges
- KERN PLS/PLS-A: force compensation
- Protective cover included in delivery
- Plug-in power supply KERN YKA-05 suitable for PLS-A and PLJ-A
- Plug-in power supply KERN YKA-04 suitable for PLS-F and PLJ-F
- To make the right choice of your KERN scale easier for you, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect scale for your application.

Additional Information

| | |
|-------------------------|--------------------------------|
| No. | KP-5396 |
| Manufacturer (Brand) | KERN |
| VGKL number | 105833620 |
| EAN | 4058072149185 |
| Transport temperature | Room temperature |
| Color | light greylgrey |
| Readability | 10 mg |
| Operation via | Membrane keypad |
| Weight | 7.9 kg |
| Reproducibility | 10 mg |
| Scales features | with switchable weighing units |
| Material weighing plate | Stainless steel |
| Weighing range MAX | 6 kg |
| Manufacturer | CORE |
| Type scale | Precision scale |

