




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN precision scale PLS, 720 g

**€1,150.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

- PLJ 2000-3A high-quality milligram balance with enormous weighing capacity up to 2100 g - ideal for large samples or heavy tare vessels
- Ergonomically optimized control panel for left- and right-handers
- KERN PLJ: internal automatic adjustment, guarantees high accuracy and makes independent of location.
- Ideal for mobile, legal-for-trade applications, such as ambulatory gold and jewelry purchasing
- KERN PLS: CAL adjustment program for setting accuracy, external test weights at extra cost
- Draft shield standard on models with weighing plate size Ø 110 mm, weighing chamber Ø×H 60×150 mm
- KERN PLS/PLJ-F: Strain gauges
- KERN PLS/PLS-A: force compensation
- Protective cover included in delivery
- Plug-in power supply KERN YKA-05 suitable for PLS-A and PLJ-A
- Plug-in power supply KERN YKA-04 suitable for PLS-F and PLJ-F
- To make the right choice of your KERN scale easier for you, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect scale for your application.

## Additional Information

No.	KP-5408
Manufacturer (Brand)	KERN
VGKL number	105831720
EAN	4058072149291
Transport temperature	Room temperature
Color	light greylgrey
Readability	1 mg
Operation via	Membrane keypad
Weight	8.1 kg
Reproducibility	1 mg
Scales features	with switchable weighing units
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Precision scale
Weighing range MAX	720 g

