




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale PNJ, 12 kg

€870.00
plus VAT &
Shipping

Product Images



Description

- pNJ only: Internal automatic calibration, guarantees high accuracy and makes independent of location. Ideal for mobile, legal-for-trade applications, such as ambulatory gold and jewelry purchasing
- pNS only: CAL adjustment program for setting accuracy using an external test weight
- High-quality tuning fork weighing system for fast display of the weight value, precise dosing and high mechanical ruggedness
- Capacity display: A rising luminous band indicates the still available weighing range
- Accurate counting: Automatic reference optimization gradually improves the average value of the part weight
- Compact dimensions, advantageous where space is limited
- Large, shock-resistant stainless steel weighing plate
- Large glass draft shield with 3 sliding doors for convenient access to the sample, for models with weighing plate size Ø 140 mm, weighing chamber W×D×H 172×171×160 mm
- Dust cover included in delivery
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model, which is made up of two quality features and, in addition to the product features, technical data and pictograms, provides you with a further decision-making aid for the perfect scale for your application.

Additional Information

| | |
|-------------------------|--------------------------------|
| No. | KP-5412 |
| Manufacturer (Brand) | KERN |
| VGKL number | 105837506 |
| EAN | 4058072149338 |
| Transport temperature | Room temperature |
| Color | light greylgrey |
| Readability | 100 mg |
| Operation via | Membrane keypad |
| Weight | 4.9 kg |
| Reproducibility | 100 mg |
| Scales features | with switchable weighing units |
| Manufacturer | CORE |
| Material weighing plate | Stainless steel |
| Type scale | Precision scale |
| Weighing range MAX | 12 kg |

