



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**KERN precision scale with battery 600 g: 0.01,
weighing plate 105 mm Ø**

**€180.00
plus VAT &
Shipping**

Product Images



Description

The scale is ready to use immediately including 9 V block batteries.

Operating time approx. 12 h.

AUTO-OFF function for battery conservation, can be switched off.

- Weighing function simplifies formulation
- particularly flat design
- safe and non-slip installation due to rubber feet
- large LCD display, digit height 15 mm
- Weighing plate made of plastic

Dimensions:

170 x 240 x 39 mm (W x D x H)

Net weight:

approx. 0.5 kg

Permissible ambient temperature range:

5°C / 35°C

Accessories: Battery 6LR61 block (Order No. 2-3027)

Power supply unit (Order No. 5-3078)

Additional Information

No.	KP-7301
Manufacturer (Brand)	KERN
EAN	4058072284183
Transport temperature	Room temperature
Color	White
Wide	240 mm
Height	39 mm
Length	170 mm
Depth	170 mm
Readability	0,0099 g
Operation via	Membrane keypad
Weight	0.5 kg
Reproducibility	0,0099 g
Scales features	ultra-flat with switchable weighing units with internal calibration with below-balance weighing system batteries/battery included
DM Weighing plate	105 mm
Energy supply	Battery Accumulator Mains
Shape weighing plate	around
Material weighing plate	Acrylonitrile butadiene styrene (ABS)
Calibration type	at the touch of a button
Weighing range MAX	600 g
Spare part	2-3027
Type scale	Precision scale
Weighing units	g oz ozt dwt

