



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Sartorius precision balance MSA12201S-000-D0 max. 12200g, 0.1 g resolution

€3,692.00
plus VAT &
Shipping

Product Images



Description

Sartorius Cubis® laboratory balance

MSA precision balance

In addition to the standard Q-Guide operator guidance system, Cubis® offers the possibility of using customized solutions with the new Q-Apps. A large number of Q-Apps are available for download. You can adapt them individually to your process flow or have them configured for your special wishes and requirements right from the start.

Cubis® Display & Control Units

MSA - The ultimate

- Top-class technology and information design
- Touchscreen with high-resolution color TFT for brilliant display of text and graphics
- Outstanding ease of use and highest display quality, especially for complex applications
- Q-Apps can be adapted to your individual workflow

Applications

Unit conversion, SQmin function for minimum weighing according to USP, automatic calibration/adjustment function isoCAL, individual labeling, density determination, statistics, calculations, averaging, formulation, weighing in percent, time-controlled functions, totalizing, DKD uncertainty of measurement, second tare memory, counting, checkweighing, alibi memory, audit trail

Leveling: automatic, manual

without draft shield

Additional Information

No.	DK-0411
Manufacturer (Brand)	neoLab
EAN	4066292504721
Transport temperature	Room temperature
Readability	0.1000 g
Operation via	Touchscreen
Shape weighing plate	angular
Weighing range MAX	12200 g
Dimensions weighing plate	208 x 208 mm
Type scale	Precision scale

