



## KERN Set EMB,

### Product Images

---



## Description

- Simple and convenient 2-button operation, thus also well suited for school and teaching use
- Weighing function simplifies formulation
- Particularly flat design
- Safe and non-slip placement due to rubber feet
- Stackable for space-saving storage
- Practical transport lock prevents overloading of the weighing system during storage or transport
- Note: The models are only supplied in sets of 5. I.e. the prices given in the table refer to a delivery of 5 pieces each. Individual delivery not possible. The indicated calibration prices refer to the calibration of a single scale
- KERN EMB 3000-1S: not in a 5-pack
- To make it easier for you to make the right choice of your KERN balance, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect balance for your application.

## Additional Information

No.	KP-4047
Manufacturer (Brand)	KERN
EAN	4058072137229
Color	hellgrau/rau
Wide	145 mm
Height	205 mm
Depth	205 mm
Operation via	Folientastatur
Weight	2.5 kg
Reproducibility	2 mg
Scales features	mit umschaltbaren Wä geeinheiten
Weighing range MAX	2 kg
Manufacturer	KERN
Type scale	Laborwaage

