



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Set EMB,

**Price on
request
plus VAT &
Shipping**

Product Images

Description

- Simple and convenient 2-button operation, thus also well suited for school and teaching use
- Weighing function simplifies formulation
- Particularly flat design
- Safe and non-slip placement due to rubber feet
- Stackable for space-saving storage
- Practical transport lock prevents overloading of the weighing system during storage or transport
- Note: The models are only supplied in sets of 5. I.e. the prices given in the table refer to a delivery of 5 pieces each. Individual delivery not possible. The indicated calibration prices refer to the calibration of a single scale
- KERN EMB 3000-IS: not in a 5-pack
- To make it easier for you to make the right choice of your KERN balance, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect balance for your application.

Additional Information

No.	KP-4050
Manufacturer (Brand)	KERN
EAN	4058072137243
Transport temperature	Room temperature
Color	light grey grey
Wide	145 mm
Height	205 mm
Depth	205 mm
Readability	0,01 g
Operation via	Membrane keypad
Weight	2.5 kg
Reproducibility	200 mg
Scales features	with switchable weighing units
Weighing range MAX	200 g
Manufacturer	CORE
Type scale	Laboratory balance

