



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN Set EMB,

**Price on  
request  
plus VAT &  
Shipping**

## Product Images

## Description

---

- Simple and convenient 2-button operation, thus also well suited for school and teaching use
- Weighing function simplifies formulation
- Particularly flat design
- Safe and non-slip placement due to rubber feet
- Stackable for space-saving storage
- Practical transport lock prevents overloading of the weighing system during storage or transport
- Note: The models are only supplied in sets of 5. I.e. the prices given in the table refer to a delivery of 5 pieces each. Individual delivery not possible. The indicated calibration prices refer to the calibration of a single scale
- KERN EMB 3000-IS: not in a 5-pack
- To make it easier for you to make the right choice of your KERN balance, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect balance for your application.

## Additional Information

---

No.	KP-4050
Manufacturer (Brand)	KERN
EAN	4058072137243
Transport temperature	Room temperature
Color	light greylgrey
Wide	145 mm
Height	205 mm
Depth	205 mm
Scales features	with switchable weighing units
Readability	0,01 g
Operation via	Membrane keypad
Weight	2.5 kg
Manufacturer	CORE
Reproducibility	200 mg
Type scale	Laboratory balance
Weighing range MAX	200 g

