



KERN Set EMB,

Product Images



Description

- Simple and convenient 2-button operation, thus also well suited for school and teaching use
- Weighing function simplifies formulation
- Particularly flat design
- Safe and non-slip placement due to rubber feet
- Stackable for space-saving storage
- Practical transport lock prevents overloading of the weighing system during storage or transport
- Note: The models are only supplied in sets of 5. I.e. the prices given in the table refer to a delivery of 5 pieces each. Individual delivery not possible. The indicated calibration prices refer to the calibration of a single scale
- KERN EMB 3000-1S: not in a 5-pack
- To make it easier for you to make the right choice of your KERN balance, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect balance for your application.

Additional Information

No.	KP-4059
Manufacturer (Brand)	KERN
EAN	4058072137298
Color	hellgraulgrau
Operation via	Folientastatur
Weight	2.5 kg
Reproducibility	200 mg
Scales features	mit umschaltbaren Wägeeinheiten
Weighing range MAX	500 g
Manufacturer	KERN
Type scale	Laborwaage

