



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN pocket scale TGC, 500 g

**€107.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- High-resolution pocket scale with particularly fast response and settling time, enabling very efficient operation
- Particularly flat design
- Cover made of impact-resistant plastic as pressure and dust protection. Can also be used as a weighing vessel
- Weighing plate stainless steel, hygienic and easy to clean
- Note: The models are only supplied in sets of 5. I.e. the prices given in the table refer to a delivery of 5 pieces each. Individual delivery not possible. The calibration prices indicated refer to the calibration of a single scale

## Additional Information

---

No.	KP-5522
Manufacturer (Brand)	KERN
EAN	4058072150303
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	100 mg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Pocket scale
Weighing range MAX	500 g

