



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN pocket scale TGD, 50 g

€280.00
plus VAT &
Shipping

Product Images



Description

- Simple and convenient 4-key operation
- Switchable from g to ct, gn, dwt, ozt, oz at the touch of a button
- Innovative touch screen: High-contrast, backlit and touch-sensitive touch display for convenient operation and comfortable reading
- Hinged lid as pressure and dust protection
- Stainless steel weighing plate, hygienic and easy to clean
- Weighing vessel as standard
- Powered by the USB cable, so no additional power supply is needed.
- Delivery in beautiful decorative packaging
- Note: The models are supplied only in sets of 5. I.e. the prices given in the table refer to a delivery of 5 pieces each. Individual delivery not possible. The calibration prices given refer to the calibration of a single scale

Additional Information

No.	KP-5524
Manufacturer (Brand)	KERN
EAN	4058072150310
Transport temperature	Room temperature
Color	light gray dark gray
Readability	1 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	50 g
Manufacturer	CORE
Type scale	Pocket scale

