



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FCB, 3 kg

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Second display on the rear of the scale PRE-TARE function for manual pre-drawing of a known container weight, useful for fill quantity checks High mobility: Thanks to battery operation/rechargeable battery operation (optional), compact, flat design and low dead weight, suitable for use at several locations (production, warehouse, shipping, etc.) RS-232 data interface as standard Dust cover included in delivery

Additional Information

| | |
|-------------------------|----------------------|
| No. | KP-4181 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072138301 |
| Transport temperature | Room temperature |
| Color | light gray dark gray |
| Readability | 100 mg |
| Operation via | Membrane keypad |
| Energy supply | Power grid |
| Material weighing plate | Stainless steel |
| Weighing range MAX | 3 kg |
| Manufacturer | CORE |
| Type scale | Bench scale |

