



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FFN, 3 kg

€325.00
plus VAT &
Shipping

Product Images



Description

Suitable for the increased hygienic requirements in the food industry, pharmaceutical industry and chemical industry
Suitable for rough industrial use Supports you in your HACCP-compliant quality system High mobility: thanks to battery operation, compact design and low own weight suitable for use in several locations such as production, kitchen, sales room, etc. Stainless steel design of housing and weighing plate. Easy to clean thanks to smooth surfaces Recessed grips on the underside for convenient transport of the scale Dust and splash protection IP65 (only in battery mode)

Additional Information

| | |
|-------------------------|----------------------|
| No. | KP-4200 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072138493 |
| Transport temperature | Room temperature |
| Color | light gray dark gray |
| Readability | 500 mg |
| Operation via | Membrane keypad |
| Energy supply | Power grid |
| Material weighing plate | Stainless steel |
| Weighing range MAX | 3 kg |
| Manufacturer | CORE |
| Type scale | Bench scale |

