



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FKB, 15 kg

**Price on
request
plus VAT &
Shipping**

Product Images

Description

- PRE-TARE function for manual pre-drawing of a known container weight, useful for fill quantity checks
- Especially practical: thanks to large weighing ranges and compact dimensions, heavy loads can be weighed with high precision even in confined spaces. Useful in determining the smallest weight differences such as consumed gas, abrasion in mechanical parts, rock samples, minerals, druses, silver, etc.
- Freely programmable weighing unit, e.g. display directly in wire length g/m, weight per unit area g/m² (paper, fabric, sheet metal, etc.) or similar.
- Protective cover included in scope of delivery

Additional Information

No.	KP-4251
Manufacturer (Brand)	KERN
EAN	4058072138844
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	1 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	15 kg
Manufacturer	CORE
Type scale	Bench scale

