



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN bench scale FKB, 16 kg

**€495.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Ideal for high-precision weighing in industrial environments thanks to the high resolution of up to 360000 points Robust die-cast aluminum housing, ensures a stable stand, protects the weighing technology and is insensitive to daily use Particularly practical: thanks to large weighing ranges and compact dimensions, heavy loads can be weighed with high precision even in the most confined spaces. Useful in determining the smallest weight differences such as consumed gas, abrasion in mechanical parts, rock samples, minerals, druses, silver, etc. PRE-TARE function for manual pre-weighing of a known container weight, useful for fill quantity checks Freely programmable weighing unit, e.g., display directly in wire length g/m, g/m2 weight per unit area (paper, cloth, sheet metal, etc.), or the like Protective work cover included in scope of delivery

## Additional Information

---

No.	KP-4253
Manufacturer (Brand)	KERN
VGKL number	105837614
EAN	4058072138868
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	100 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	16 kg
Manufacturer	CORE
Type scale	Bench scale

