



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FKB, 65 kg

**Price on
request
plus VAT &
Shipping**

Product Images

Description

Ideal for high-precision weighing in industrial environments thanks to the high resolution of up to 360000 points Robust die-cast aluminum housing, ensures a stable stand, protects the weighing technology and is insensitive to daily use Particularly practical: thanks to large weighing ranges and compact dimensions, heavy loads can be weighed with high precision even in the most confined spaces. Useful in determining the smallest weight differences such as consumed gas, abrasion in mechanical parts, rock samples, minerals, druses, silver, etc. PRE-TARE function for manual pre-weighing of a known container weight, useful for fill quantity checks Freely programmable weighing unit, e.g., display directly in wire length g/m, g/m² weight per unit area (paper, cloth, sheet metal, etc.), or the like Protective work cover included in scope of delivery

Additional Information

No.	KP-4258
Manufacturer (Brand)	KERN
VGKL number	105837618
EAN	4058072138912
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	500 mg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Bench scale
Weighing range MAX	65 kg

