



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FOB, 500 g

€125.00
plus VAT &
Shipping

Product Images



Description

Stainless steel design of housing and weighing plate, thus rustproof and easy to clean thanks to smooth surfaces High mobility: thanks to compact, flat design suitable for use in several locations (kitchen, sales room, market stall, food laboratory, etc.) Particularly flat design Safe and non-slip installation thanks to rubber feet IP65: protected against dust and splash water (only in battery operation) Ready to start: Battery 9 V block as standard, operating time up to 24 h Protective work cover included in scope of delivery

Additional Information

No.	KP-4287
Manufacturer (Brand)	KERN
EAN	4058072139193
Transport temperature	Room temperature
Color	light gray dark gray
Readability	100 mg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Bench scale
Weighing range MAX	500 g

