



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN bench scale FOB, 15 kg

**€335.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

NEW: Innovative weighing with tolerance range (checkweighing): The color of the display changes depending on the weighing result (too light/ok/too heavy) and thus supports portioning, dosing and sorting Stainless steel design of housing and weighing plate. Easy to clean thanks to smooth surfaces Increased moisture protection thanks to waterproof silicone sealing of load cell, electronics and soldering points High mobility: thanks to compact, flat design suitable for use in several locations (kitchen, sales room, market stall, food laboratory, etc.) Safe and non-slip installation thanks to rubber feet Control of the weighing range already occupied at the touch of a button (net/gross weight) Supports you in your HACCP-compliant quality system Protective work hood included in delivery Application examples: Weighing ingredients and spices in the catering industry Portioning dough, meat, fish, poultry, salad plates in canteens, etc. Filling of cans etc. Checking of prepackages

## Additional Information

---

No.	KP-4288
Manufacturer (Brand)	KERN
EAN	4058072139209
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	5 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	15 kg
Manufacturer	CORE
Type scale	Bench scale

