



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FOB, 8 kg | 15 kg

€295.00
plus VAT &
Shipping

Product Images



Description

- Stainless steel housing and weighing plate, thus rust-proof and easy to clean thanks to smooth surfaces
- Supports you in your HACCP-compliant quality system
- High mobility: thanks to compact, flat design suitable for use in several locations (kitchen, sales room, market stall, food laboratory, etc.)
- Safe and slip-free installation due to rubber feet
- Dust and splash water protection IP67 (only in battery operation)
- Ready to start: Batteries 4×1.5 V AA as standard, operating time up to 48 h
- Protective work cover included

Additional Information

No.	KP-4289
Manufacturer (Brand)	KERN
EAN	4058072139216
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	1 g
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Bench scale
Weighing range MAX	15 kg

