



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FOB, 1,500 kg

€240.00
plus VAT &
Shipping

Product Images



Description

NEW: Innovative weighing with tolerance range (checkweighing): The color of the display changes depending on the weighing result (too light/ok/too heavy) and thus supports portioning, dosing and sorting Stainless steel design of housing and weighing plate. Easy to clean thanks to smooth surfaces Increased moisture protection thanks to waterproof silicone sealing of load cell, electronics and soldering points High mobility: thanks to compact, flat design suitable for use in several locations (kitchen, sales room, market stall, food laboratory, etc.) Safe and non-slip installation thanks to rubber feet Control of the weighing range already occupied at the touch of a button (net/gross weight) Supports you in your HACCP-compliant quality system Protective work hood included in deliveryApplication examples: Weighing ingredients and spices in the catering industry Portioning dough, meat, fish, poultry, salad plates in canteens, etc. Filling of cans etc. Checking of prepackages

Additional Information

No.	KP-4290
Manufacturer (Brand)	KERN
EAN	4058072139223
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	500 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	1 kg
Manufacturer	CORE
Type scale	Bench scale

