



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FOB, 5 kg

€125.00
plus VAT &
Shipping

Product Images



Description

- Stainless steel design of housing and weighing plate, thus rustproof and easy to clean thanks to smooth surfaces
- High mobility: thanks to compact, flat design suitable for use in several locations (kitchen, sales room, market stall, food laboratory, etc.)
- Particularly flat design
- Safe and slip-free installation due to rubber feet
- IP65: Protected against dust and splash water (only in battery operation)
- Ready to start: Battery 9 V block as standard, operating time up to 24 h
- Protective work cover included in delivery

Additional Information

| | |
|-------------------------|----------------------|
| No. | KP-4299 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072139292 |
| Transport temperature | Room temperature |
| Color | light gray dark gray |
| Readability | 1 g |
| Operation via | Membrane keypad |
| Energy supply | Power grid |
| Manufacturer | CORE |
| Material weighing plate | Stainless steel |
| Type scale | Bench scale |
| Weighing range MAX | 5 kg |

