



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FOB, 5 kg | 7,500 kg

€295.00
plus VAT &
Shipping

Product Images



Description

- Stainless steel housing and weighing plate, thus rust-proof and easy to clean thanks to smooth surfaces
- Supports you in your HACCP-compliant quality system
- High mobility: thanks to compact, flat design suitable for use in several locations (kitchen, sales room, market stall, food laboratory, etc.)
- Safe and slip-free installation due to rubber feet
- Dust and splash water protection IP67 (only in battery operation)
- Ready to start: Batteries 4×1.5 V AA as standard, operating time up to 48 h
- Protective work cover included

Additional Information

No.	KP-4301
Manufacturer (Brand)	KERN
EAN	4058072139315
Transport temperature	Room temperature
Color	light gray dark gray
Readability	500 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	0 kg
Manufacturer	CORE
Type scale	Bench scale

