



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FXN, 30 kg

€460.00
plus VAT &
Shipping

Product Images



Description

- Suitable for the increased hygienic requirements in the food industry
- Supports you in your HACCP compliant quality system
- Suitable for harsh industrial use
- IP68: very high protection class against dust and water. Therefore very robust and durable. Ideal for harsh environmental conditions. Even submersion of the scale in water possible
- Stainless steel design of housing and weighing plate, thus rust-proof and easy to clean thanks to smooth surfaces
- Very fast display: stable weighing values within approx. 2 sec
- High mobility: thanks to battery operation, compact design and low weight, suitable for use at several locations (production, kitchen, sales room, etc.)
- Recessed grips on the underside for convenient transport of the scale
- Checkweighing: a visual and acoustic signal supports portioning, dosing, and sorting

Additional Information

No.	KP-4314
Manufacturer (Brand)	KERN
EAN	4058072139438
Transport temperature	Room temperature
Color	light gray dark gray
Readability	5 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	30 kg
Manufacturer	CORE
Type scale	Bench scale

