



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale FXN, 3 kg

€460.00
plus VAT &
Shipping

Product Images



Description

- Suitable for the increased hygienic requirements in the food industry
- Supports you in your HACCP compliant quality system
- Suitable for harsh industrial use
- IP68: very high protection class against dust and water. Therefore very robust and durable. Ideal for harsh environmental conditions. Even submersion of the scale in water possible
- Stainless steel design of housing and weighing plate, thus rust-proof and easy to clean thanks to smooth surfaces
- Very fast display: stable weighing values within approx. 2 sec
- High mobility: thanks to battery operation, compact design and low weight, suitable for use at several locations (production, kitchen, sales room, etc.)
- Recessed grips on the underside for convenient transport of the scale
- Checkweighing: an optical and acoustic signal supports portioning, dosing and sorting

Additional Information

No.	KP-4316
Manufacturer (Brand)	KERN
EAN	4058072139452
Transport temperature	Room temperature
Color	light gray dark gray
Readability	500 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	3 kg
Manufacturer	CORE
Type scale	Bench scale

