



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN Bench Scale GAB, 6 kg | 15 kg

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- Compact dimensions, advantageous where space is limited
- High mobility: thanks to battery operation (optional) and compact, flat design, suitable for use at several locations (production, warehouse, shipping, etc.)
- Weighing with tolerance range (checkweighing): an optical signal supports portioning, dosing and sorting
- Totalizing of weight values
- Protective work cover included in the scope of delivery

Examples of use:

- Defection scales for pharmacies
- Hand mixing of tea, coffee, chocolates
- Portioning of dough, meat, fish, poultry, salad plates in canteens etc.
- Weighing crops ambulatory by kilogram price
- Checkweighing in supermarkets
- High precision industrial applications
- **Note:** Official calibration is mandatory for commercial use

## Additional Information

---

No.	KP-4321
Manufacturer (Brand)	KERN
EAN	4058072139490
Transport temperature	Room temperature
Color	light gray dark gray
Readability	2 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	15 kg
Manufacturer	CORE
Type scale	Bench scale

