



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN Bench Scale GAB, 6 kg

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

Compact dimensions, advantageous where space is limited High mobility: thanks to battery operation (optional) and compact, flat design, suitable for use at several locations (production, warehouse, shipping, etc.) Weighing with tolerance range (checkweighing):an optical signal supports portioning, dosing and sorting Totalizing weight values Protective work cover included in scope of deliveryApplication examples: Defection scales for pharmacies Manual mixing of tea, coffee, chocolates Portioning of dough, meat, fish, poultry, salad plates in canteens, etc. Ambulatory weighing of crops by kilogram price Checkweighing in supermarkets High-precision industrial applicationsNote: Official calibration is mandatory in commercial transactions

## Additional Information

No.	KP-4325
Manufacturer (Brand)	KERN
EAN	4058072139537
Transport temperature	Room temperature
Color	light gray dark gray
Readability	50 mg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Bench scale
Weighing range MAX	6 kg

