



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Bench Scale GAB, 3 kg | 6 kg

Price on
request
plus VAT &
Shipping

Product Images



Description

- Compact dimensions, advantageous where space is limited
- High mobility: thanks to battery operation (optional) and compact, flat design, suitable for use at several locations (production, warehouse, shipping, etc.)
- Weighing with tolerance range (checkweighing): an optical signal supports portioning, dosing and sorting
- Totalizing of weight values
- Protective work cover included in the scope of delivery

Examples of use:

- Defection scales for pharmacies
- Hand mixing of tea, coffee, chocolates
- Portioning of dough, meat, fish, poultry, salad plates in canteens etc.
- Weighing crops ambulatory by kilogram price
- Checkweighing in supermarkets
- High precision industrial applications
- **Note:** Official calibration is mandatory for commercial use

Additional Information

No.	KP-4326
Manufacturer (Brand)	KERN
EAN	4058072139544
Transport temperature	Room temperature
Color	light gray dark gray
Readability	1 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	6 kg
Manufacturer	CORE
Type scale	Bench scale

