



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN bench scale WTB 1.5 kg

€285.00
plus VAT &
Shipping

Product Images



Description

- Suitable for the increased hygienic requirements in the food industry
- Supports you in your HACCP compliant quality system
- Dust and splash water protection IP65 (only in battery mode)
- Very fast display: stable weighing values within approx. 2 sec
- Second display on the back of the scale
- Compact dimensions, advantageous where space is limited
- Innovative weighing with tolerance range (checkweighing): The color of the display changes depending on the weighing result (too light/ok/too heavy) and thus supports portioning, dosing and sorting
- High mobility: thanks to battery operation (optional) and compact, flat design suitable for use at several locations
- Recessed grips on the underside for convenient transport of the scale
- Protective work cover included in the scope of delivery

Additional Information

No.	KP-5584
Manufacturer (Brand)	KERN
EAN	4058072150853
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	200 mg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Bench scale
Weighing range MAX	1 kg

