



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN bench scale WTB 3000 g

**€285.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Suitable for the increased hygienic requirements in the food industry
- Supports you in your HACCP compliant quality system
- Dust and splash water protection IP65 (only in battery mode)
- Very fast display: stable weighing values within approx. 2 sec
- Second display on the back of the scale
- Compact dimensions, advantageous where space is limited
- Innovative weighing with tolerance range (checkweighing): The color of the display changes depending on the weighing result (too light/ok/too heavy) and thus supports portioning, dosing and sorting
- High mobility: thanks to battery operation (optional) and compact, flat design suitable for use at several locations
- Recessed grips on the underside for convenient transport of the scale
- Protective work cover included in the scope of delivery

## Additional Information

---

No.	KP-5589
Manufacturer (Brand)	KERN
EAN	4058072150907
Transport temperature	Room temperature
Color	light gray dark gray
Readability	500 mg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Bench scale
Weighing range MAX	3 kg

