



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN counting system Max 300 kg: $d=0.00001$ kg

€1,170.00
plus VAT &
Shipping

Product Images



Description

With this highly accurate counting system KERN CCS a wide range of single counting scales can be replaced cost-effectively and efficiently. Reference scale KERN CFS: This professional counting scale KERN CFS, which can also be used individually, meets even the highest demands on accuracy, weighing range and piece volumes by connecting a high-load weighing bridge Programmable via keypad: - desired reference piece count - known reference weight Three displays for weight display, reference weight, total piece count Memory (PLU) for 100 articles with additional text and tare weight e.g. of a container for 100 articles with additional text, reference weight and tare weight, e.g. of a container Acoustic fill-to-target: target piece count or target weight programmable. When the target value is reached, a signal sounds Accurate counting: The automatic reference optimization gradually improves the average value of the part weight Quantity scale KERN KFP: The quantity piece counting takes place with high precision on the weighing platform (= weighing bridge) KERN CCS. So even the smallest counting parts can be counted in largest volumes Stainless steel platform, base steel painted Aluminum single point load cell (1 x 3000 e), dust and splash water protection IP65 Delivery time for models with >= 600kg: 5 working days (freight forwarding) Protective work cover included in delivery

Additional Information

No.	KP-2181
Manufacturer (Brand)	KERN
EAN	4058072129835
Transport temperature	Room temperature
Readability	0,001 g
Energy supply	Power grid
Weighing range MAX	300 kg
Manufacturer	CORE
Type scale	Counting scale

