



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN counting scale CFS, 300 g

€510.00
plus VAT &
Shipping

Product Images



Description

- Memory (PLU) for 100 articles with additional text, reference weight and tare weight of e.g. a container
- Accurate counting: Automatic reference optimization gradually improves the average value of the part weight
- Checkweighing: a visual and acoustic signal supports portioning, dosing and sorting
- Programmable via keypad
 - - desired reference piece count
 - - known reference weight
- Three displays for weight indication reference weight, total piece count
- Fill-to-target function
- Programmable target piece count or target weight, e.g. for checkweighing. Reaching the target value is indicated by a signal
- PRE-TARE function for manual preweighing of a known container weight, useful for fill quantity checks
- Second scale interface for setting up a high-resolution counting system, e.g. with KFP V20 weighbridges
- Dust cover included in the scope of delivery

Additional Information

No.	KP-0095
Manufacturer (Brand)	KERN
EAN	4058072124571
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	1 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	300 g
Manufacturer	CORE
Type scale	Counting scale

