



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN counting scale CFS, 6 kg

€340.00
plus VAT &
Shipping

Product Images



Description

Memory (PLU) for 100 items with additional text, reference weight and tare weight of e.g. a container Accurate counting: Automatic reference optimization gradually improves the average value of the part weight Checkweighing: an optical and acoustic signal supports portioning, dosing and sorting Programmable via keypad: - desired reference piece count - known reference weight Three displays for weight indication Reference weight, total piece count Fill-to-target function Target piece count or target weight programmable, e.g. for checkweighing. Reaching the target value is indicated by a signal PRE-TARE function for manual preweighing of a known container weight, useful for fill quantity checks Second scale interface for setting up a high-resolution counting system, e.g. with KFP V20 weighbridges Dust cover included in delivery

Additional Information

| | |
|-------------------------|----------------------|
| No. | KP-0099 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072124618 |
| Transport temperature | Room temperature |
| Color | light gray dark gray |
| Readability | 100 mg |
| Operation via | Membrane keypad |
| Energy supply | Power grid |
| Material weighing plate | Stainless steel |
| Weighing range MAX | 6 kg |
| Manufacturer | CORE |
| Type scale | Counting scale |

