



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN counting scale CKE, 3,600 kg

€425.00
plus VAT &
Shipping

Product Images



Description

- Self-explanatory graphic-supported control panel, counting process immediately comprehensible even without operating instructions.
- - no training time = saves costs
- - ideal for the inexperienced user
- - visualized sequence avoids operating errors
- The 4 steps are carried out from left to right:
- Place empty container on the weighing plate and tare by pressing a key (TARE)
- Fill the container with the reference sample quantity (5, 10 or 20 pieces)
- Confirm the selected reference sample quantity by pressing the key (5, 10 or 20)
- Fill the counted quantity into the container. The piece count is shown directly on the display
- Accurate counting: Automatic reference optimization gradually improves the average value of the part weight
- Two scales in one: switch from counting mode to weighing mode at the touch of a button
- Battery operation possible, only for models with weighing plate size W×D 340×240 mm, batteries 6×1.5 V Size C not standard, operating time up to 40 h
- Protective work cover included in delivery

Additional Information

No.	KP-0713
Manufacturer (Brand)	KERN
EAN	4058072126391
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	10 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	3 kg
Manufacturer	CORE
Type scale	Counting scale

