



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN CPB counting scale, 6 kg | 15 kg

€290.00
plus VAT &
Shipping

Product Images



Description

- Accurate counting: Automatic reference optimization gradually improves the average value of the part weight
- Programmable via keypad:
 - - desired reference piece count
 - - known reference weight
- Three displays for weight indication (verifiable), reference weight, total piece count
- Counting total memory: for totalizing equal counting parts in total piece and total weight
- Acoustic fill-to-target: target piece count or target weight programmable, e.g. for checkweighing. A signal sounds when the target value is reached
- PRE-TARE function for manual preweighing of a known container weight, useful for fill quantity checks
- High mobility: thanks to battery operation (optional), compact, flat design and low weight, suitable for use at several locations (production, warehouse, shipping, etc.)
- Two scales in one: switch from counting mode to weighing mode at the touch of a button
- Protective work cover included in delivery
- Note: Official calibration is mandatory for commercial use

Additional Information

No.	KP-0318
Manufacturer (Brand)	KERN
EAN	4058072125189
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	2 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	15 kg
Manufacturer	CORE
Type scale	Counting scale

