



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN CPB counting scale, 30 kg

€290.00
plus VAT &
Shipping

Product Images



Description

Accurate counting: Automatic reference optimization gradually improves the average value of the part weight
Programmable via keypad: – desired reference piece count – known reference weight
Three displays for weight indication (verifiable), reference weight, total piece count
Counting total memory: for adding up equal counting parts in total piece and total weight
Acoustic fill-to-target: target piece count or target weight programmable, e.g. for checkweighing. When the target value is reached, a signal sounds
PRE-TARE function for manual pre-drawing of a known container weight, useful for fill quantity checks
High mobility: thanks to battery operation (optional), compact, flat design and low dead weight, suitable for use at several locations (production, warehouse, shipping, etc.)
Two scales in one: switch from counting mode to weighing mode at the touch of a button
Protective work cover included in delivery
Note: Official verification is mandatory for commercial use

Additional Information

No.	KP-7290
Manufacturer (Brand)	KERN
old neoLab article no.	KP-0319
EAN	4058072284084
Transport temperature	Room temperature
Color	light gray dark gray
Readability	500 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	30 kg
Manufacturer	CORE
Type scale	Counting scale

