



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN CXB counting scale, 30 kg

€350.00
plus VAT &
Shipping

Product Images



Description

- Accurate counting: Automatic reference optimization gradually improves the average value of the part weight

Programmable via keypad:

- desired reference piece count
- known reference weight
- Three displays for weight indication (verifiable), reference weight, total piece count
- Acoustic fill-to-target: target piece count or target weight programmable, e.g. for checkweighing. A signal sounds when the target value is reached
- Counting total memory: for totalizing equal counting parts in total piece and total weight
- Energy management: backlight off after 5 sec
- PRE-TARE function for manual pre-weighing of a known container weight, this weight can be stored in memory
- Two scales in one: switch from counting mode to weighing mode at the touch of a button
- Dust cover included in delivery
- **Note:** Official verification is mandatory for commercial use

Additional Information

No.	KP-3958
Manufacturer (Brand)	KERN
VGKL number	105833002
EAN	4058072136369
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	2 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	30 kg
Manufacturer	CORE
Type scale	Counting scale

