



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN CXB counting scale, 3 kg

**€350.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Accurate counting: Automatic reference optimization gradually improves the average value of the part weight Programmable via keypad: – desired reference piece count – known reference weight Three displays for weight indication (verifiable), reference weight, total piece count Acoustic fill-to-target: target piece count or target weight programmable, e.g. for checkweighing. When the target value is reached, a signal sounds Counting total memory: for adding up equal counted parts in total piece and total weight Energy management: backlighting off after 5 sec PRE-TARE function for manual preweighing of a known container weight, this can be stored Two scales in one: switchover from counting mode to weighing mode at the touch of a key Protective work cover included in delivery Note: Official verification is mandatory in commercial transactions

## Additional Information

---

No.	KP-3959
Manufacturer (Brand)	KERN
VGKL number	105830602
EAN	4058072136376
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	200 mg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Counting scale
Weighing range MAX	3 kg

