



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Binder water supply set, for treatment or demineralization of tap water

€1,201.00
plus VAT &
Shipping



Description

External water supply set for the treatment or full desalination of tap water.
Consisting of fresh and waste water canisters of 20 liters each, wiring and pump.
To be hung on the back of the unit.

Additional Information

No.	BB-0285
Manufacturer (Brand)	Binder
old neoLab article no.	7-5812
EAN	4058072513566
Transport temperature	Room temperature
Display (reading, setting)	graphic display
Operation via	Knobs/Buttons
Application product water	General laboratory applications (e.g. autoclaving, media preparation, etc.)
Withdrawal possibilities product water type II	Canister
Model High purity/clean water system	Binder Constant climatic chambers
Quality product water	Deionized water (quality depends on feed water)
Feed water quality	Drinking water quality
Volume water reservoir	20 l

