



Julabo PURA 10 water bath

Product Images



Description

PURA water baths stand for innovative, simple and safe working in the laboratory. All models (except PURA 4) are supplied with an integrated drain screw for quick and clean emptying of the baths. PURA water baths can be used flexibly and are easy to transport thanks to ergonomic recessed grips on the sides. Non-slip rubber feet ensure a firm footing and more safety in the laboratory. Cooling coil: Accessories Bath vessel: Stainless steel Bath cover: Accessories Medium bath opening (cm/WxD/BD): 22 x 27 / 14 Product features: - Easy operation - Luminous display - Splash-proof power switch - Built-in boil-dry protection - Integrated drain plug (except PURA 4) - No disturbing elements in the bath - Adjustable shelf supplied as standard - Integrated timer - Non-slip rubber feet - Ergonomic recessed grips Performance data: For temperature control tasks near or below ambient temperature: use cooling coil or JULABO immersion cooler. Temperature stability with bath cover (accessory)

Additional Information

No.	JN-0711
Manufacturer (Brand)	Julabo
VGKL number	103590510
EAN	4058072848552
Net weight	7.00000 kg
Country of origin	Deutschland
Customs tariff number	84198998

