



Julabo PURA 4 water bath

Product Images



Description

PURA water baths stand for innovative, simple and safe working in the laboratory. All models (except PURA 4) are supplied with an integrated drain screw for quick and clean emptying of the baths. PURA water baths can be used flexibly and are easy to transport thanks to ergonomic recessed grips on the sides. Non-slip rubber feet ensure a firm footing and more safety in the laboratory. Cooling coil: Accessories Bath vessel: Stainless steel Bath cover: Accessories Medium bath opening (cm/WxD/BD): 12 x 27 / 14 Product features: - Easy to operate - Luminous display - Splash-proof power switch - Built-in boil-dry protection - Integrated drain plug (except PURA 4) - No disturbing elements in the bath - Shelf supplied as standard - Integrated timer - Non-slip rubber feet - Ergonomic recessed grips Performance data: Temperature constancy with bath cover (accessory)

Additional Information

No.	JN-0710
Manufacturer (Brand)	Julabo
VGKL number	103590504
EAN	4058072848545
Net weight	4.90000 kg
Country of origin	Deutschland
Customs tariff number	84198998

