



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



sunlab® stainless steel water bath, 20 liters

€1,115.00
plus VAT &
Shipping

Product Images



Description

The **Sunlab® stainless steel water baths** are the ideal solution for heating samples and media. The temperature range can be set from RT +5 °C to 99 °C. Temperature setting is quick and easy via the digital control panel. After use, the water bath can be easily emptied via a bottom outlet (D-8811 and D-8812).

- with digital temperature display
- with limit value alarm

Additional Information

No.	D-8812
Manufacturer (Brand)	sunlab
EAN	4058072227852
Transport temperature	Room temperature
Temp MAX	95 °C
Wide	400 mm
Height	370 mm
Length	572 mm
Depth	572 mm
Operation via	Membrane keypad
Weight	18 kg
Voltage	230 V
Water baths and thermostats properties	with overheating protection with closure/lid
Energy supply	Power grid
Heating capacity	1200 W
Material bath	Stainless steel
Bath depth	140 mm
Bath opening	270 x 450 mm
Temp. accuracy	+/- 0.3 °C
Volume water bath	20 l
Material inside	Stainless steel
Temp MIN above ambient temperature	5 °C
Type temperature control	PID

